

Technical Data

CROS Pure™



Data Sheet

www.bestsound-technology.com

SIEMENS

CROS Pure

Style Specific Features


Ingress Protection Rating	IP67
Nanocoated housing	●
RoHS	Yes

Battery Specifications

Battery type	312
Battery voltage	1.3 V
Battery current drain	1.4 mA
Battery life (cell zinc air)	~90 h

Technical Data

IRIL IEC 118-13:2011 (bystander)	
800-960 MHz	<-5 dB SPL
1400-2000 MHz	<-10 dB SPL
Field strength magnetic link	
EN300330 short range device	-24.5 dB μ A/m @ 10 m
RS 210, issue 8	-17.5 dB μ V/m @30m
Sending frequency	3.28 MHz
Transmission range	18 cm

 **WARNING**

Choking hazard posed by small parts.

- ▶ This instrument is not intended for the fitting of infants, small children and persons of mental incapacity.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice. The required features should therefore be specified in each individual case at the time of conclusion of the respective contract.

© 11.2015 AR, Siemens AG
Order No. 02794-99T1-7600

Sivantos GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany
Phone +49 9131 308 0